

ComfiLife

Adjustable Ankle Sleeve

THANK YOU!

We appreciate your trust in our products, we are ecstatic that you chose ComfiLife Adjustable Ankle Sleeve over many other possibilities out there. You will soon realize that you have made the best decision!

Now sit back and relax, we've got your back!

Not Satisfied?

Contact us

Satisfaction Guaranteed

At ComfiLife we take pride in the quality of our products and our customer service. Your purchase comes with a 100% satisfaction and replacement guarantee.

If you are not 100% satisfied for any reason or if there are any issues with your ankle sleeve please contact us and we will send you a FREE REPLACEMENT.

For any questions please contact us and our lightning speed customer service will take care of you.

Please Read Before Use

- Your ankle sleeve fits left or right foot. It is normal to feel some tightness around the ankle during the first few uses. The fit should be snug and not too tight, it should not restrain you from doing your daily activities.
- Have rest periods between uses to let your skin breathe to avoid skin irritation.
- The ankle sleeve straight from the package may have a slight odor. This is normal, air it out for a few hours and wash before use. Constant use may stretch the material of your ankle brace, washing it would restore the firmness of the elastic fabric.
- Always consult with your physician before using an ankle support brace.

Other Best Sellers From ComfiLife

ComfiLife Ankle Support Wrap

ComfiLife Ankle Sleeve

ComfiLife Knee Sleeve

ComfiLife Hot & Cold Pack

ComfiLife Knee Ice Pack



ComfiLife COMFORT ANYWHERE

How to use?

Left Foot



Place your **left** foot at the center of the sleeve marked "**left**" aligning with the arrow.



Bottom of your heel should be partially covered.



Wrap the thicker straps around your ankle.



Fasten the thicker strap on each other.



Pull the v-shaped elastic bands up.



Stretch the bands to the desired level of compression.



Fasten the bands to the right and left sides of your ankle.



Make sure the fit is snug but not too tight.

Right Foot



Place your **right** foot at the center of the sleeve marked "**right**" aligning with the arrow.



Bottom of your heel should be partially covered.



Wrap the thicker straps around your ankle.



Fasten the thicker strap on each other.



Pull the v-shaped elastic bands up.



Stretch the bands to the desired level of compression.



Fasten the bands to the right and left sides of your ankle.



Make sure the fit is snug but not too tight.

Care Instructions

• Hand wash in cold water. Wash gently, do not scrub. Do not use fabric softener, bleach or any harsh chemicals. Air dry and make sure your wrap is completely dry before putting it on.

WARNING: DO NOT PUT YOUR WRAP IN THE DRYER OR ON A HEATER. DO NOT USE HARSH DETERGENTS OR BLEACH PRODUCTS.



Hand Wash Cold



Do Not Bleach



Do Not Dry Clean



Air Dry



Do Not Iron

Visit: www.ComfiLife.com

within 14 days to register your product for full warranty