

# THANK YOU!

We appreciate your trust in our products, we are ecstatic that you chose ComfiLife Knee Sleeve over many other possibilities out there. You will soon realize that you have made the best decision!

Now sit back and relax, we've got your back!

## Not Satisfied?

Contact us

### Satisfaction Guaranteed

At ComfiLife we take pride in the quality of our products and our customer service. Your purchase comes with a 100% satisfaction and replacement guarantee.

If you are not 100% satisfied for any reason or if there are any issues with your knee sleeve please contact us and we will send you a FREE REPLACEMENT.

For any questions please contact us and our lightning speed customer service will take care of you.

### Please Read Before Use

- Your knee sleeve fits left or right knee. It is normal to feel some tightness around the knee during the first few uses. The fit should be snug and not too tight, it should not restrain you from doing your daily activities.
- Have rest periods between uses to let your skin breathe to avoid skin irritation.
- The knee sleeve straight from the package may have a slight odor. This is normal, air it out for a few hours and wash before use. Constant use may stretch the material of your knee sleeve, washing it would restore the firmness of the elastic fabric.
- Always consult with your physician before using a knee sleeve.

### Other Best Sellers From ComfiLife

ComfiLife Ankle Support Wrap



ComfiLife Knee Brace



ComfiLife Patella Brace



ComfiLife Hot & Cold Pack



ComfiLife Knee Ice Pack



## How to use?

The ComfiLife Knee Sleeve's Compression Fit Technology provides optimal knee support and promotes muscle recovery.

Suitable for sports activities with its breathable and flexible material.

Wear your knee sleeve such that the horizontal lines that form an oblong shape sits on your knee cap. Also make sure that the ComfiLife Logo is not upside down while you wear it.



## Care Instructions

- Hand wash in cold water. Wash gently, do not scrub. Do not use fabric softener, bleach or any harsh chemicals. Air dry and make sure your sleeve is completely dry before putting it on.

**WARNING: DO NOT PUT YOUR SLEEVE IN THE DRYER OR ON A HEATER. DO NOT USE HARSH DETERGENTS OR BLEACH PRODUCTS.**



Do Not Dry Clean



Do Not Bleach



Hand Wash Cold



Air Dry



Do Not Iron

Visit: [www.ComfiLife.com](http://www.ComfiLife.com)

**within 14 days to register your product for full warranty**